

Utility

Category: Pattern

Approvals

- ✓ UX, Harry Ardon
- ✓ SG, Jayson Elliot
- ✓ ADA, Mark Stimpson
- ✓ FED, Quentin Franke

Status

- ✓ In Web Style Guide
- ✓ In (Development or Production)
- ✓ Shared AEM Component
- ✗ Not in Mobile App

UTILITY

1. Symbols

Symbol Library

Help |  Print

Help

Icons

Do this



Don't do this




Arrow





Alert




Annotations


 Type something

Type something 

 Type something


Type something 

Screens



2. Examples

Example 1


Elena ▾

My Health Medical Record Message Center Appointments **Pharmacy** Coverage & Costs Health & Wellness
🔍

Pharmacy

Request and track prescription refills for yourself and those you care for. 🛒 Fill order (0)

Elena Martinez ▾
Help | 🖨️ Print

Ready for fill (4)

Rx

Levothyroxine

Last filled 4/7/2018

Status: 2 refills remaining

[Drug information](#) | [Add to order](#)

Rx

Adalimumab

Last filled 4/3/2018

Status: 5 refills remaining

[Drug Information](#) | [Add to order](#)

Rx

Spiriva Handihaler

Last filled 3/21/2018

Status: 3 refills remaining

[Drug Information](#) | [Add to order](#)

Rx

Rosuvastatin

Last filled 3/17/2018

Status: 2 refills remaining

[Drug Information](#) | [Add to order](#)

More ways to get care

📊

Rx estimate

[Get a drug cost estimate ↗](#)

📄

Refill by Rx number

[Find prescriptions for yourself or someone else](#)

📍

Pharmacy locator

[Look for a pharmacy in your area](#)

📄

Transfer Rx

[Transfer prescriptions to a Kaiser Permanente pharmacy](#)

💬

Chat with a pharmacist

Weekdays, 8 a.m. - 5 p.m.

[Start a chat now ↗](#)

More resources

[Drug encyclopedia](#)

[Drug formulary](#)

All other prescriptions

Rx

Levothyroxine

Last filled 4/7/2018

Status: 2 refills remaining

[Drug Information](#)

Rx

Adalimumab

Last filled 4/3/2018

Status: 5 refills remaining

[Drug Information](#)

Rx

Spiriva Handihaler

Last filled 3/21/2018

Status: 3 refills remaining

[Drug Information](#)

Rx

Rosuvastatin

Last filled 3/17/2018

Status: 2 refills remaining

[Drug Information](#)

1-5 of 19 items

<
1
2
3
4
>

Get Care

- [Advice](#)
- [Routine Care](#)
- [Urgent Care](#)
- [Emergency Care](#)
- [Find Doctors & Locations](#)
- [Behavioral Health](#)
- [Health Classes](#)
- [When Traveling](#)

Our Organization

- [About KP](#)
- [News & Views](#)
- [Commitment to the Community](#)
- [Diversity & Inclusion](#)
- [Awards & Accreditations](#)
- [Annual Report](#)
- [Careers](#)
- [Media Inquiries](#)

Member Support

- [New Member Welcome](#)
- [Forms & Publications](#)
- [Member Assistance](#)
- [Member Services](#)

Visit Our Other Sites

- [Individual & Family Plans](#)
- [Medicaid/Medi-Cal](#)
- [Medicare](#)
- [Affordable Care Act](#)
- [For Businesses](#)
- [Broker Support](#)

Language

- [Español](#)
- [🌐 Other Languages](#)







Accessibility
Privacy
Terms & Conditions
Rights & Responsibilities
Site Policies
Disaster Planning
Technical Information

Site Map
Contact Site Manager

Copyright © 2016 Kaiser Foundation Health Plan, Inc.



You will need the free [Adobe Acrobat Reader](#) to read PDFs.




UTILITY


2. Examples

Example 2


 **Medical Record** 

Bobbie Sue Zorkendörfer 

[< Medical Record](#) [Help](#)

Test Results 

Click on a test name to view details.
Please note that not all of your test results may be listed.

Sort: Name A-Z 

LABEL COL
[Lorem Ipsum](#)

LABEL COL
[Lorem Ipsum](#)

LABEL COL
[Shock and/or unconsciousness](#)

LABEL COL
[Lorem Ipsum](#)

LABEL COL
[Nausea and/or vomiting](#)

LABEL COL
[Lorem Ipsum](#)



LABEL COL
[Skin rash and/or hives](#)

LABEL COL
[Lorem Ipsum](#)

LABEL COL
[Sulfa](#)

LABEL COL
[Skin rash and/or hives](#)

1-5 of 25 items






 of 10 pages 


Certain content delivered by MyChart [®], licensed from Epic Systems Corporation, © 1999 to 2014, patents pending.


You may not see all the tests ordered by your doctor on this page. Some test results may still be pending at this time, and some test results cannot be released electronically because of regulatory restrictions.


If you have questions about these test results, please e-mail your doctor or call 404-365-0966.


Follow us

Get Care 


Our Organization 

Member Support 

Visit Our Other Sites 

[Español](#)

© 2016 Kaiser Foundation Health Plan, Inc.



UTILITY

3. Usage Guidelines

3.1 What is utility?

The Utility Component contains the links “Help” and “Print”. “Help” links to the help section and “Print” allows a user to print the contents of a page.

There are certain screens within each ACAA center where the Utility Component is hidden or modified. Please refer to each individual center’s documentation for more details.


Print functionality is custom to the Center where the print link manifests— see individual Center/Page specs for descriptions of Print functionality on case-by-case basis.

For general Print functionality specs, see “Print” section of UX Style Guide.


UTILITY

4. Design Guidelines

Desktop



Tablet




A Consistent UI for both desktop and tablet



At 600 and below, only "Help" exists. "Print" is omitted at this breakpoint.

(320 shown as example)

Desktop




Tablet



1px divider Hex: E9EBEE

600 and below (320 shown as example)




Functionality

Clicking "Help" opens the help page in the same window
Clicking "Print" opens the print page in a new window

Styling

Desktop and tablet





Mobile




Alignment & spacing

Desktop and tablet



Mobile



A Utility component is right aligned with last column on the grid.

Next UI element is spaced 40px away from utility item. Could be text, a button, etc.

UTILITY

5. Interaction Guidelines

5.1 Interaction Guidelines

There are certain screens within each ACAA center where the Utility Component is hidden or modified. Please refer to each individual center's documentation for more details.

In most other cases, the Utility Component follows the rules below.

Print functionality is custom to the Center where the print link manifests— see individual Center/Page specs for descriptions of Print functionality on case-by-case basis. For general Print functionality specs, see “Print” section of UX Style Guide.

Help and Print exist at the 768 breakpoint and above.


Clicking “Help” opens the help page in the same window

Clicking “Print” opens the print page in a new window


5. Interaction Guidelines

5.2 Print output variations - Desktop / Tablet


Variation 1



Header and Nav removed





Header and Nav removed




Text should be default to the original design (hex color), the user should have the option to print in black from their printer settings.

5. Interaction Guidelines

5.3 Print Flow



*PRINT OCCURS




UTILITY

6. Development Guidelines

6.1 KP Pattern Library


<http://cskpcloudxn1713.cloud.kp.org/patterns-library/print-utility.html>

6.2 Desktop Redlines



Gotham Styling 6
Book
KP Blue

6.3 Mobile Redlines



Gotham Styling 7
Book
KP Blue

UTILITY

7. ADA Considerations

7.1 ADA

This Pattern follows standard coding guidelines for accessibility.

UTILITY

8. Version

11.19.2019

Pattern documented

UTILITY

9. Depricated